

The crowning in level of detail

AFFINIS
PERFECT IMPRESSIONS



Still one of the fastest A-silicones

AFFINIS fast, with a maximum total setting time of 2 min 10 sec guarantees optimum efficiency and saves time.



Precise moulding of every detail even under the most difficult conditions

The highly active surfactant system is characterized by excellent surface affinity in a moist environment.




Confidence during impression taking

Self-contouring consistency prevents air bubbles and distortions.

















Easy monitoring for success

AFFINIS PRECIOUS wash materials with impressive gold and silver pigmentation for excellent readability and detail reproduction.







AFFINIS Family

Impression Material

REF	Product	Consistency	Max. processing time	Intraoral setting time	Total setting time	Handmix	System 50/75	System 360
6501	light body		01:00	02:00	03:00		50ml ●	
6601	light body fast		00:50	01:20	02:10		50ml ●	
6511	regular body		01:00	02:00	03:00		50ml ●	
6611	regular body fast		00:50	01:20	02:10		50ml ●	
6775	PRECIOUS light body		01:00	02:20	03:20		50ml ●	
6776	PRECIOUS regular body		01:00	02:20	03:20		50ml ●	
6547	MonoBody		01:00	02:30	03:30		75ml ●	
6482			01:00	02:30	03:30			380ml Refill ●
6520	heavy body		01:00	02:00	03:00		75ml ●	
6487			01:00	02:00	03:00			380ml Refill ●
6620	heavy body fast		00:50	01:20	02:10		75ml ●	
60019554	Black Edition (heavy body)		01:00	02:20	03:20		75ml ●	
60019776			01:00	02:20	03:20			380ml Refill ●
6535	putty super soft		01:05	02:00	03:05	600ml		
6530	putty soft		01:10	02:00	03:10	600ml		
6531	putty soft fast		01:05	01:50	02:55	600ml		
6472	putty		01:00		03:30			380ml Refill ●

Accessories

REF	Product		pcs
60019663	Universal Mixing Tip (12x ●)		40
60019694			100
6552	Mixing Tip Turquoise (8x ●)		40
6555	Oral Tips Yellow		100
6162			Dynamic Mixing Tips Yellow
6163	Fixation Ring (1x ●)		1

© COLTENE – www.coltene.com



Imprint / Publisher

Coltene/Whaledent AG
Feldwiesenstrasse 20
CH-9450 Altstätten
+41 71 757 5300
info.ch@coltene.com

Sales Contact

Coltene/Whaledent Ltd.
The President Suite
Kendal House, Victoria Way
Burgess Hill, West Sussex
RH15 9NF/UK
+44 1444 235486
info.uk@coltene.com

Service Center

service@coltene.com
00 800 11 0000 12 (UIFN)
Mon.-Fri. 8:00 to 17:00 MESZ

 COLTENE